

Azide activation kit

100-250 µg

INSTRUCTIONS

Last revised May 2020

Instructions for product no:

L1-AZ1-025

GlyCLICK® Azide activation of up to 250 µg IgG

Product description

GlyCLICK Azide activation is for site specific azide activation of up to 250 µg of antibody. The conserved N-linked glycosylation site on the CH2 domain of each heavy chain of the Fc region is used by GlyCLICK for selective antibody modification. The azide activation kit introduces GalNAz site-specifically on Fc for subsequent click reaction.

Immobilized GlycINATOR® removes all Fc N-glycans, including high-mannose, hybrid-type and bisected glycans to the inner GlcNAc. The subsequent azide activation at the GlcNAc enables a following click reaction using e. g. copper free strain promoted azide-alkyne cycloaddition (SPAAC) for selective attachment of a label of choice. The modification primes the antibody for conjugation of the desired molecule at activated sites on the Fc regions for incorporation of two labels per antibody (DOL=2), see Figure 1.

The modification procedure is performed by combining two enzymatic steps at the Fc domain of the IgG. All steps are performed under physiological conditions, thus maintaining the quality of the antibody. The site-specific modification on the Fc domain preserves the affinity of the antigen binding sites.




Figure 1. The modification is performed in two steps:

1. Immobilized GlycINATOR hydrolyzes the N-glycans on the Fc-part of the IgG to the inner GlcNAc.
2. Azide attachment on the GlcNAc using GalT(Y289L)* and UDP-GalNAz*.

*GalT(Y289L) and UDP-GalNAz are components of SiteClick™ and are provided under an intellectual property license from Life Technologies Corporation. The trademark SiteClick™ is the property of Life Technologies Corporation.

Customer Service USA Phone: 617-444-8421 | info@genovis.com | Genovis AB, 245 First Street, Cambridge, MA 02142, USA
Customer Service Europe Phone: 0046 (0)46 10 12 30 | info@genovis.com | Genovis AB, Box 790, SE-220 07 Lund, Sweden
Support support@genovis.com **Order** order@genovis.com | Fax: 0046 (0)46 12 80 20

www.genovis.com

Content and storage

GlyCLICK Azide activation kit contains enzymes, reagents and material to azide activate 100-250 µg antibody.

GlyCLICK Azide activation kit is shipped cold and components should be stored at 4 °C to 8 °C upon arrival. **Note: Do NOT freeze!**

Table 1. Content and storage temperatures of GlyCLICK components.


Name	Amount	Store at
Antibody concentrator (small), 50K	1 piece	4 °C to 25 °C
Collection tube for small concentrator	2 pieces	4 °C to 25 °C
Desalting Spin column, 40K	1 piece	4 °C to 8 °C
Immobilized GlycINATOR [®] , spin column	1 piece	4 °C to 8 °C
UDP-GalNAz (0.22 mg)	1 vial solid	4 °C to 8 °C Protect from light
20× TBS pH 7.4 (0.5 M)	1.8 mL	4 °C to 8 °C
Buffer additive	30 µL	4 °C to 8 °C Protect from light
β-1,4-galactosyltransferase (Y289L) (GalT)	25 µL	4 °C to 8 °C Protect from light
Antibody concentrator (large), 50K	1 piece	4 °C to 25 °C

Before you begin, briefly centrifuge tubes. Always wear suitable laboratory protective clothing and gloves when handling these reagents.

Do not freeze any of the components of the kit!

GlyCLICK Azide activation kit is for R&D use only.

Overview of the protocol for antibody azide activation using GlyCLICK



Equipment required:

- Centrifuge with fixed angle rotor that can accommodate 1.5-2 ml centrifuge tubes.
- Centrifuge with swinging bucket rotor that can accommodate 17 mm × 100 mm centrifuge tubes.
- Incubator or water bath for 25 °C and 30 °C.
- End-over-end mixer

Additional Materials required

- 100 to 250 µg of antibody, max volume of 200 µl in a non-phosphate or Tris-based buffer, free of carrier proteins and/or azide. 20x TBS, a desalting spin column (40K) for buffer exchange and a small concentrator (50K) is provided for convenience. For adjusting the antibody solution, please follow "Guidance for concentration and buffer exchange" below.
- Centrifuge tubes: 1.5-2 ml and 15 ml.
- ddH₂O. **Note: if a chelating agent will be used as label it is important to use metal free water (trace analysis grade) throughout the protocol.**
If labeling is performed with DIBO-DFO, or other conjugation reagent with chelator, the antibody must not be in contact with glass or metal.

Sodium azide must be avoided throughout the protocol!

Guidance for concentration and buffer exchange

The antibody concentration step **is required if**:

- The volume of the antibody is more than 200 μl .

If the sample volume is 100-200 μl but needs a buffer exchange (if it contains phosphate or azide), follow the instruction in section "Buffer exchange with Desalting column, 0.5 ml". It is advisable to start with more antibody than 250 μg if concentration or buffer exchange of the sample is needed prior to "Step 1. Modification of the carbohydrate on antibody Fc domain, deglycosylation".

Concentration step

1. Add 500 μl of ddH₂O to the small antibody concentrator and cap the device as shown in Figure 2.
2. Centrifuge at 5000 \times g for 6 minutes. Make sure that **the cap strap and one membrane panel of the concentrator face the center of the rotor** (Fig. 2).
3. Discard the flow-through.
4. Add the antibody solution to the small antibody concentrator.
5. Centrifuge at 5000 \times g for 2-6 minutes. Make sure that **the cap strap and one membrane panel of the concentrator face the center of the rotor** (Fig. 2).
6. Discard the flow-through.

Note: If the antibody volume in the concentrator is more than 200 μl , centrifuge for an additional 2 minutes at 5000 \times g, or until the appropriate volume is achieved.

7. Invert the small antibody concentrator into the collection tube as shown in Figure 2.
8. Centrifuge at 1000 \times g for 3 minutes to collect the concentrated antibody. After collection, the amount of concentrated Ab should be approximately 150-200 μl in the collection tube.




Figure 2. Antibody concentration step

If you need to buffer exchange the sample, follow the steps below.

Buffer exchange with Desalting Spin column, 0.5 ml

This step **is required if**:

- The antibody is in a phosphate-based buffer (e.g. PBS), and/or
- The antibody is in a buffer containing azide

1. Prepare 10 ml 1x TBS buffer by adding 500 μl 20x TBS to 9.5 ml ddH₂O in a 15 ml tube. Vortex briefly to mix.
2. Break off the bottom closure of the Desalting Spin column. Loosen the lid (**do not** remove the lid).
3. Place the column in a collection tube (1.5-2 ml) and centrifuge at 1500 \times g for 1 min to remove the storage solution.
4. Discard the flow-through and place the column in the collection tube.

5. Add 300 μ l 1x TBS buffer on top of the resin. Centrifuge the column at 1500 \times g for 1 min and discard the flow-through.
6. Repeat step 5 **two more times**. Last spin for 2 minutes.
7. Blot the bottom of the column to remove excess liquid. Place the column in a new tube (1.5-2 ml).
8. Apply the antibody solution on top of the resin (100-200 μ l).
9. Centrifuge at 1500 \times g for 2 min and retain the flow-through containing the antibody in 1x TBS buffer.

Detailed protocol for azide activation of 100 - 250 μ g of antibody

Step 1. Modification of the carbohydrate on Antibody Fc domain, deglycosylation

The antibody solution should be in non-phosphate based buffer with no azide. Max 250 μ g in 200 μ l.

Time Required: 15 minutes hands-on, 30 minutes hands-off

Materials from kit: 1x TBS buffer (prepared from 20x TBS),
Spin column with Immobilized GlycINATOR[®]

- The lid and the cap of the spin column are used during the incubation
 - Before the centrifugations, remove the bottom cap and slightly open the lid
 - Let the Immobilized GlycINATOR column equilibrate to room temperature before use.
- 1.1 Prepare 10 ml 1x TBS buffer by adding 500 μ l 20x TBS to 9.5 ml ddH₂O in a 15 ml tube. Vortex briefly to mix.
 - 1.2 Break off the bottom plastic cap of the GlycINATOR column (save the cap) and slightly open the lid. Place the column in a 1.5-2 ml collection tube.
 - 1.3 Centrifuge the column at 200 \times g for 1 min to remove the storage solution.
 - 1.4 Discard the flow-through.
 - 1.5 Place the column in the collection tube.
 - 1.6 Add 300 μ l 1x TBS buffer to the top of the resin. Centrifuge the column at 200 \times g for 1 minute and discard the flow-through.
 - 1.7 Repeat the steps in 1.6 **two more times**.
 - 1.8 Re-insert the bottom cap at the bottom of the spin column.
 - 1.9 Immediately add the antibody solution (100- 200 μ l) to the column. Re-seal the column with the lid.
 - 1.10 Beware; fully suspend the media manually and make sure it is flowing in the column.
 - 1.11 Incubate the column by end-over-end mixing at room temperature for 30 minutes.
 - 1.12 Remove the bottom cap and place the column in a clean micro centrifuge tube (1.5-2 ml). Loosen the lid.
 - 1.13 Centrifuge the column at 1000 \times g for 1 minute to collect the deglycosylated antibody sample.

For maximum recovery of the sample:

- 1.14 Attach the bottom cap. Add 50 μ l 1x TBS and seal the column with the lid.
- 1.15 Invert the column a couple of times.
- 1.16 Remove the bottom cap and place the column in a clean micro centrifuge tube (1.5-2 ml). Loosen the lid.
- 1.17 Centrifuge the tube at 1000 \times g for 1 minute to collect the deglycosylated antibody sample.
- 1.18 Repeat steps 1.14 to 1.17 once more.
- 1.19 Pool the collected deglycosylated antibody samples.

Step 2. Azide attachment

Time required: 5 minutes hands-on, followed by overnight incubation

Materials from kit: 1× TBS buffer (prepared from 20× TBS),
UDP-GalNAz,
GalT enzyme,
Buffer additive

- 2.1 Prepare the azide modification solution by adding the following components to the tube containing UDP- GalNAz.

Add to the UDP-GalNAz tube:

- 30 µL of buffer additive
- Deglycosylated Ab solution (from step 1.13 and 1.19) and 1x TBS buffer to a total of 375 µL.

Mix the solution by carefully pipetting up and down.

- 2.2 Add the GalT enzyme, 25 µL. The final reaction volume should be 400 µL. Mix the solution by carefully pipetting up and down. Wrap the tube cap with Parafilm® or similar.
- 2.3 Incubate overnight, protected from light at 30 °C.

Step 3. Purification and concentration of azide activated Antibody

Time required: 1 hour

Materials from kit: 1× TBS buffer (prepared from 20x TBS),
Large antibody concentrator

- This step will remove excess of UDP-GalNAz

Wash antibody concentrator

- 3.1 Prepare 10 ml of 1× TBS by adding 500 µL of 20× TBS to 9.5 ml of ddH₂O in a 15 ml tube. Vortex briefly to mix.
- 3.2 Remove the conical collection tube from the large antibody concentrator as shown in Figure 3.
- 3.3 Add 2 ml of 1× TBS to the large antibody concentrator and centrifuge at 1200 × g for 10 minutes. **Make sure that one membrane panel of the concentrator faces the center of the rotor.**
- 3.4 Discard the flow-through.

Purify the antibody

- 3.5 Add 1.6 ml of 1× TBS and 400 µL of the azide activated antibody from Step 2.3 to the large antibody concentrator (Fig. 3).
- 3.6 Centrifuge at 1200 × g for 6 minutes. **Make sure that one membrane panel of the concentrator faces the center of the rotor.**
- 3.7 Discard the flow-through.
- 3.8 Add 1× TBS to a total of 2 ml to the large antibody concentrator.

- 3.9 Centrifuge at $1200 \times g$ for 10 minutes. **Make sure that one membrane panel of the concentrator faces the center of the rotor.**
- 3.10 Discard flow-through.
- 3.11 Repeat steps 3.8 -3.10 **two more times.**

Note: If the antibody volume in the concentrator is more than $\sim 200 \mu\text{l}$, the volume in the concentrator can be reduced by additional centrifugation e.g. for an additional 5 minutes at $1200 \times g$ or until the appropriate volume is achieved.




Figure 3. Use of large concentrator

Antibody collection

- 3.12 Invert the antibody concentrator into the conical collection tube as shown in Figure 4.
- 3.13 Centrifuge at $1000 \times g$ for 3 minutes to collect the concentrated azide activated antibody.
- 3.14 Transfer the azide activated antibody from the conical collection tube to a 1.5 ml centrifuge tube.
- 3.15 If Nanodrop is available, determine protein concentration.
- 3.16 At this stage, the azide activated antibody can be stored at $2-8 \text{ }^\circ\text{C}$ for conjugation of label at a later time.




Figure 4. Collection of purified and concentrated azide activated antibody

References

1. Sjögren, J. et al., 2013. EndoS2 is a unique and conserved enzyme of serotype M49 group A Streptococcus that hydrolyses N-linked glycans on IgG and α 1-acid glycoprotein. *The Biochemical Journal*, 455(1), pp.107–118.
2. Ramakrishnan, B. & Qasba, P.K., 2002. Structure-based design of beta 1,4-galactosyltransferase I (beta 4Gal-T1) with equally efficient N-acetylgalactosaminyltransferase activity: point mutation broadens beta 4Gal-T1 donor specificity. *J Biol Chem*, 277(23), pp.20833–20839.

Legal and Disclaimers

All rights reserved. Genovis products are covered by one or more patents, trademarks and/or copyrights owned or controlled by Genovis AB. For more information about commercial rights, please contact Genovis team at info@genovis.com. Genovis products are intended for research use only. They are not intended to be used for therapeutic or diagnostic purposes in humans or animals. All goods and services are sold subject to Genovis' General Terms and Conditions of Sale.

SiteClick™

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.

The trademark SiteClick™ is the property of Life Technologies Corporation.

©2020 Genovis AB.