
FabRICATOR®

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www.genovis.com

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SmartEnzymes™

GENOVIS

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FabRICATOR® 2000 units (A0-FR1-020)

Digestion of up to 2.0 mg IgG

FabRICATOR® 5000 units (A0-FR1-050)

Digestion of up to 5.0 mg IgG

FabRICATOR® 5x5000 units (A0-FR1-250)

Digestion of up to 5 × 5.0 mg IgG

FabRICATOR® Validation kit (A0-FR4-060)

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FabRICATOR®

Limited Use Label License: Research Use Only

The purchase of the **IdeS** enzyme from *Streptococcus pyogenes* (sold under the trade name FabRICATOR®) conveys to the purchaser the limited, non-transferable right to use the purchased amount of **IdeS** only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel.

Purchaser agrees to be bound by the following terms and restrictions:

1) A right is granted purchaser only for internal research purposes using **IdeS** for digesting an IgG and is not for use in commercial services of any kind, including, without limitation, reporting the result of purchaser's activities for a fee or other form of consideration.

2) **IdeS** will not be made available by purchaser to any third parties in any form, separately or in combination, for any monetary or other consideration or at no charge, except that **IdeS** may be made available to third parties who agree to be bound by all the terms and restrictions of this right for purposes of evaluation only.

3) **IdeS** and the digested IgG will not be used *in vivo* in humans.

4) Purchaser will not make commercial use of the **IdeS** unless it first

secures a Sublicense Agreement from Genovis AB for such commercial use.

All rights reserved. For information on obtaining additional rights, please contact: licensing@genovis.com

Genovis AB holds an exclusive world wide license to all patents derived from international publication WO03051914, including granted US Patent No US 7,666,582 regarding IdeS from *Streptococcus pyogenes* for biotechnical industrial applications which are neither therapeutic nor diagnostic, other than the following exception which is included within the license: digesting IgG *in vitro* in clinical samples for diagnostic purposes.

This license is defined in the intellectual property license agreement between Hansa Biopharma AB and Genovis AB.

The trademark FabRICATOR® is the property of Genovis AB.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

All goods and services are sold subject to Genovis' General Terms and Conditions of Sale.

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© Genovis AB
Genovis Inc.