



FucosEXO®

Immobilized

STORE AT

+4-8°C



FOR RESEARCH USE ONLY

Instructions for Use

FucosEXO® Immobilized

Microspin 5 × 0.5 mg (G1-FM6-025)

Process 5 × 0.5 mg glycoprotein

FucosEXO® Immobilized

Microspin 10 × 0.5 mg (G1-FM6-050)

Process 10 × 0.5 mg glycoprotein

DOWNLOAD INSTRUCTIONS FOR USE



www.genovis.com/ifu-G1-FM6

Immobilized Enzyme for Hydrolysis of α 1-2,3,4 Fucose in Spin Columns

FucosEXO is a mix of α -Fucosidases for efficient hydrolysis of α 1-2, α 1-3 and α 1-4-linked fucose residues on native N- and O-glycosylated proteins or free oligosaccharides. The FucosEXO Immobilized spin columns contain two fucosidases covalently coupled to agarose beads, for defucosylation of glycoproteins or free oligosaccharides without contaminating the final preparation with enzyme. FucosEXO hydrolyzes fucose on glycoproteins under native conditions and display a high activity in pH 6.0-8.0.

The enzymes in FucosEXO are derived from *Akkermansia muciniphila* and *Streptococcus oralis*, and expressed in *E. coli*. Both enzymes contain a His-tag and the molecular weights of the enzymes are 87 kDa and 64 kDa.

CONTENT AND STORAGE

The FucosEXO Immobilized columns contain sufficient material to defucosylate 0.5 mg glycoprotein per column. The resin is supplied in 20% ethanol with no preservatives added.

FucosEXO Immobilized is shipped cold and should be stored at +4-8°C upon arrival. Do not freeze the product!

FucosEXO Immobilized is for R&D use only.

QUALITY CONTROL

FucosEXO Immobilized is tested to meet the specifications and lot-to-lot consistency.

FucosEXO Immobilized is tested for absence of microbial contamination with blood agar plates, Sabouraud dextrose agar plates and fluid thioglycollate medium.

YOU MIGHT ALSO BE INTERESTED IN

FucosEXO® Lyophilized

Lyophilized enzyme mix for hydrolysis of α 1-2,3,4 fucose

SialEXO™

Hydrolysis of sialic acids

GalactEXO®

Hydrolysis of β 1-3,4 galactose

GaINAcEXO®

Hydrolysis of α -linked GaINAcS

Preparations

Important Information

- Use lids and bottom caps during the incubation.
- Before centrifugation, remove the bottom cap and loosen the lid (do not remove the lid).

Additional Materials Required

- Reaction buffer: 20mM Tris pH 6.8.^{1,2}
- Microcentrifuge tubes (1.5-2 ml).

1. FucosEXO Immobilized displays high activity in buffers with pH 6.0-8.0, and over a wide range of ionic strengths (0-500mM NaCl). Optimization might be required if a buffer other than the recommended reaction buffer is used.
2. If the glycoprotein sticks to the resin, a buffer with NaCl can be used in the reaction and/or in the wash steps.

Hydrolysis of α 1-2,3,4 Fucose in Spin Columns

Sample Preparation

Prepare the glycoprotein in 100-300 μ l reaction buffer per column. Max amount of glycoprotein is 0.5 mg per column.

1. Equilibration

- 1.1 Break off the bottom cap of the FucosEXO Immobilized column (save the cap) and place the column in a microcentrifuge tube. Loosen the lid.
- 1.2 Centrifuge at 200 \times g for 1 min to remove the storage solution. Discard the flow-through.
- 1.3 Equilibrate the column by adding 300 μ l reaction buffer and centrifuge at 200 \times g for 1 min. Discard the flow-through.
- 1.4 Perform step 1.3 two additional times.
- 1.5 Insert the bottom cap.

2. Enzymatic Reaction

- 2.1 Add the glycoprotein in a volume of 100-300 μ l reaction buffer. Max 0.5 mg glycoprotein per column.
- 2.2 Seal the column with the lid.
- 2.3 Fully suspend the media, mix by inversion and make sure there is a flow in the column.
- 2.4 Incubate the column with end-over-end mixing at room temperature for 60 min³.

3. Collection of Processed Material

- 3.1 Remove the bottom cap and place the column in a new microcentrifuge tube. Loosen the lid.
- 3.2 Centrifuge at 1000 \times g for 1 min to collect the processed material.

4. For Maximum Recovery of the Sample

- 4.1 Insert the bottom cap.
- 4.2 Add 100 μ l reaction buffer.²
- 4.3 Seal the column with the lid and invert it a couple of times.
- 4.4 Remove the bottom cap and place the column in a new microcentrifuge tube. Loosen the lid.
- 4.5 Centrifuge at 1000 \times g for 1 min to collect the processed material.
- 4.6 Repeat steps 4.1-4.5.
- 4.7 Pool the collected fractions, including the sample from step 3.2.

3. Longer incubation times may be required depending on the glycoprotein.

USA & Canada

Genovis Inc.

10919 Technology Place Suite C, San Diego, CA 92127, USA

Phone: 1-855-782-0084 (toll free)

Fax: 1-858-524-3006

EMEA & Asia

Genovis AB

Box 4, SE-24421 Kävlinge, Sweden

Phone: +46 46 10 12 30

Fax: +46 46 12 80 20

support@genovis.com

www.genovis.com



All rights reserved. Genovis products may be covered by one or more patents, trademarks and copyrights owned or controlled by Genovis AB. For more information about commercial rights, please contact the Genovis team at licensing@genovis.com.

Genovis products are intended for research use only. They are not intended to be used for therapeutic or diagnostic purposes in humans or animals.

All goods and services are sold subject to Genovis' General Terms and Conditions of Sale.

© Genovis AB