



SialEXO™

Lyophilized

STORE AT

-20°C



FOR RESEARCH USE ONLY

Instructions for Use

SialEXO™ Lyophilized 2000 units (G1-SM1-020)
Process 2 mg glycoprotein

DOWNLOAD INSTRUCTIONS FOR USE



www.genovis.com/ifu-G1-SM1

Lyophilized Enzyme Mix for Hydrolysis of Sialic Acids

The SialEXO products are sialidases for efficient desialylation of N- and O-glycosylated proteins. SialEXO Lyophilized consists of two sialidases, each having a unique activity, and they work simultaneously. One enzyme has a broad activity for α 2-3, α 2-6 and α 2-8-linked sialic acids, and the other one quickly hydrolyzes sialic acids with α 2-3-linkages. SialEXO Lyophilized removes all sialic acids on native proteins within two hours.

SialEXO hydrolyzes glycoproteins under native conditions and displays activity at pH 6.5-9.0.

The enzymes in SialEXO Lyophilized are derived from *Akkermansia muciniphila* and expressed in *E. coli*. Both enzymes contain a His-tag and the molecular weights of the enzymes are 43 kDa and 66 kDa.

UNIT DEFINITION

One unit SialEXO Lyophilized hydrolyzes sialic acids from $\geq 90\%$ of 1 μ g glycoprotein (fetuin) when incubated in 20 mM Tris pH 6.8 at 37°C for 2 hours.

CONTENT AND STORAGE

SialEXO Lyophilized is supplied lyophilized in TBS, pH 7.6, with no preservatives added.

SialEXO Lyophilized is shipped at ambient temperature, and should be stored at -20°C upon arrival.

After reconstitution, SialEXO Lyophilized is stable for at least 1 month at +4-8°C.

SialEXO is for R&D use only.

QUALITY CONTROL

SialEXO Lyophilized is tested to meet the specifications and lot-to-lot consistency.

SialEXO Lyophilized is tested for absence of microbial contamination with blood agar plates, Sabouraud dextrose agar plates and fluid thioglycollate medium.

YOU MIGHT ALSO BE INTERESTED IN

SialEXO™ Immobilized

Immobilized enzyme mix for hydrolysis of sialic acids in spin columns

OpeRATOR™

O-glycan-specific protein digestion

GalactEXO®

Hydrolysis of β 1-3,4 galactose

GalNAcEXO®

Hydrolysis of α -linked GalNAcs

FucosEXO®

Hydrolysis of α 1-2,3,4 fucose

OmniGLYZOR®

Hydrolysis of N- and mucin-type O-glycans

Preparations

Additional Materials Required

- Reaction buffer: 20mM Tris, pH 6.8.¹

1. SialEXO displays high activity in buffers at pH 6.5-9.0.

Hydrolysis of Sialic Acids

Sample Preparation

Prepare the glycoprotein in the reaction buffer. The final protein concentration in the digestion reaction should be 0.1-5 mg/ml.

1. Prepare SialEXO

1.1 Reconstitute SialEXO in 50 μ l ddH₂O to a concentration of 40 units/ μ l.

2. Add SialEXO

2.1 Add 1 unit SialEXO / 1 μ g glycoprotein or to 1 pmol oligosaccharides.²

3. Enzymatic Reaction

3.1 Incubate for 30 min-2 h at 37°C.³

2. A higher enzyme concentration may increase digestion efficiency of individual glycoproteins or oligosaccharides. This requires optimization.
3. Longer incubation times may be required depending on the glycoprotein.

USA & Canada

Genovis Inc.

10919 Technology Place Suite C, San Diego, CA 92127, USA

Phone: 1-855-782-0084 (toll free)

Fax: 1-858-524-3006

EMEA & Asia

Genovis AB

Box 4, SE-24421 Kävlinge, Sweden

Phone: +46 46 10 12 30

Fax: +46 46 12 80 20

support@genovis.com

www.genovis.com



All rights reserved. Genovis products may be covered by one or more patents, trademarks and copyrights owned or controlled by Genovis AB. For more information about commercial rights, please contact the Genovis team at licensing@genovis.com.

Genovis products are intended for research use only. They are not intended to be used for therapeutic or diagnostic purposes in humans or animals.

All goods and services are sold subject to Genovis' General Terms and Conditions of Sale.

© Genovis AB