

## INSTRUCTIONS

Version 15.1.1

Instructions for product no		
A0-IZ6-010	2 columns	Deglycosylation of up to 2x0.5mg IgG
A0-IZ6-025	5 columns	Deglycosylation of up to 5x0.5mg IgG
A0-IZ6-050	10 columns	Deglycosylation of up to 10x0.5mg IgG

### Product Description

IgGZERO® is an endoglycosidase with a very high specificity for IgG molecules of all species and subclasses.

deGlycIT™ MicroSpin contains IgGZERO® covalently coupled to agarose beads for deglycosylation of IgG. IgG is incubated with the IgGZERO® agarose beads for 15 min, deglycosylated IgG is then collected by a 1 minute centrifugation step. Since IgGZERO® is immobilized on agarose beads there is no need for extensive purification to remove the enzyme.

### Content and storage

Spin columns containing IgGZERO® covalently coupled to agarose beads.

deGlycIT™ MicroSpin are supplied in 20% EtOH and no preservatives are added.

One micro spin column contains sufficient IgGZERO® coupled agarose beads to deglycosylate 0.5 mg IgG.

deGlycIT™ MicroSpin is shipped on ice. deGlycIT™ MicroSpin should be stored at +4-8°C upon arrival.

deGlycIT™ MicroSpin is for R&D use only.

### Quality Control deGlycIT™

deGlycIT™ is tested to ensure lot-to-lot consistency.

deGlycIT™ is tested for absence of microbial contamination with blood agar plates, Sabaraud dextrose agar plates and fluid thioglycolate medium.

### Additional Materials Required

- ✓ Cleavage buffer: 10mM sodium phosphate, 150mM NaCl, pH 7.4.
- ✓ Collection tubes: Micro centrifuge tubes.

### Method

- ✓ Make sure your antibody is in cleavage buffer (See Additional Material Required above).
- ✓ Break off the bottom plastic cap of the spin column. Remember to save the bottom cap!
- ✓ Lids and bottom caps are used during the incubation.
- ✓ Before centrifugation remove the bottom cap and slightly open the lid ~90° counter clockwise.

1. Break off the bottom cap of the spin column and slightly open the screw cap lid ~90° counter clockwise.
2. Centrifuge the column at 200xg for 1min to remove storage solution.
3. Equilibrate the column with 300µl cleavage buffer (See Additional Material Required above).

4. Centrifuge the column at 200×g for 1min.
5. Repeat step 3 and 4 two times.
6. Re-insert the bottom cap into the bottom of the spin column.
7. Immediately add 100µl IgG at a maximal concentration of 5mg/ml in cleavage buffer.
8. Re-seal the column with the lid.
9. Take care to fully suspend the media manually and make sure it is flowing in the column.
10. Incubate the column by end-over-end mixing at room temperature for 15 min.
11. Twist open the lid and remove the bottom cap.
12. Place the column in a clean micro centrifuge tube (not included).
13. Centrifuge the column at 1000×g for 1min to elute the sample.

For maximum recovery of your sample:

14. Place the column in a clean micro centrifuge tube (not included).
15. Add 100µL cleavage buffer.
16. Centrifuge the column at 1000×g for 1min to elute your sample.
17. Repeat step 14-16 one more time.

©2015 Genovis AB. All rights reserved. Aspects of IgGZERO® technology are encompassed by pending patent applications in the name of Genovis AB.

The trademark IgGZERO is the property of Genovis AB.

**For research use only.** Not intended for any animal or human therapeutic diagnostic use.