

Issued 2017-02-23

Product Name	GlycINATOR®
Product Number	A0-GL1-020
Product composition	2000 Units of GlycINATOR® endoglycosidase enzyme lyophilized in Tris buffer saline (TBS).
Manufacturer	Genovis AB, Sweden
Appearance	White to light yellow powder
Source	Cloned from <i>Streptococcus pyogenes</i> and overexpressed in <i>E.coli</i>
Bioburden	Absence of microbial contamination*
Enzyme activity	One unit is defined as the amount of enzyme required to remove \geq 95% of N-linked Fc glycans from 1 μ g of human IgG1 in 10mM phosphate, 150mM NaCl pH 7.4 in 37°C for 30 minutes, as monitored by SDS-PAGE
Physical purity	Purified to \geq 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.
Shipping conditions	Ambient
Storage temp	-20°C

*Bioburden Tested on Sabouraud Dextrose Agar, Fluid Thioglycollate Medium and Blood Agar.

Telephone: +46 (0)46 10 12 30
Telefax: +46 (0)46 12 80 20
E-mail: info@genovis.com
Internet: www.genovis.com

Visiting address:
Medicon Village, Scheelevägen 2
SE-223 63 Lund
Sweden

Postal address:
Box 790
SE-223 07 Lund
Sweden

VAT number: SE556574534501

Telephone: +46 (0)46 10 12 30
Telefax: +46 (0)46 12 80 20
E-mail: info@genovis.com
Internet: www.genovis.com

Visiting address:
Medicon Village, Scheelevägen 2
SE-223 63 Lund
Sweden

Postal address:
Box 790
SE-223 07 Lund
Sweden

VAT number: SE556574534501