

Product Name	FragIT™ MicroSpin
Product Number	A0-FR6-010
Product composition	2x0.1 ml* settled resin of FabRICATOR® (IdeS) immobilized on highly cross-linked agarose. Agarose resins are suspended in 20% ethanol/water (V/V)
Manufacturer	Genovis AB, Sweden
Storage temp	4-8°C
Appearance	White resin
Source	FabRICATOR® (IdeS) enzyme is cloned from <i>Streptococcus pyogenes</i> and overexpressed in <i>E.coli</i> . Enzyme is conjugated to highly cross-linked agarose beads.
Sterility	Absence of bacterial growth**
Stability	1 year at 4-8°C

\*Immobilized enzyme activity      0.1 ml of settled resin fragments ≥95% of 0.5 mg of human IgG1 in 10 mM Phosphate, 150 mM NaCl, pH 7.4 for 15 min at 25°C as monitored by SDS-PAGE.

\*\*Sterile tests      Tested on Fluid Thioglycollate Medium and Blood Agar and Sabouraud Dextrose Agar.