



Last revised January 2020

Product Name	GlycINATOR® LE
Product Number	A0-GL8-020
Product Composition	2000 units of GlycINATOR endoglycosidase enzyme lyophilized in Tris buffer saline (TBS).
Manufacturer	Genovis AB, Sweden.
Appearance	White to light yellow powder.
Source	GlycINATOR (EndoS2) enzyme is cloned from <i>Streptococcus pyogenes</i> and expressed in <i>E. coli</i> .
Bioburden*	Absence of microbial contamination.
Enzyme Activity	One unit is defined as the amount of enzyme required to remove $\geq 95\%$ of N-linked Fc glycans from 1 μg of human IgG1 in 10mM phosphate, 150mM NaCl pH 7.4 in 30 minutes at 37°C as monitored by SDS-PAGE.
Physical Purity	Purified to $\geq 95\%$ homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.
Endotoxin level**	< 0.2 EU/vial
Shipping Conditions	Shipped at ambient temperature.
Storage Temperature	-20°C

* Tested on Sabouraud Dextrose Agar, Fluid Thioglycollate Medium and Blood Agar.

** Tested according to current EU pharm 2.6.14 method D at certified laboratory.