



Last revised February 2020

Product Name	Anti-FabRICATOR® Antibody
Product Number	A3-AF1-010
Product Description	Goat polyclonal antibody to FabRICATOR
Purity	Protein G purified
Isotype	IgG
Concentration	4 mg/ml
Immunogen	FabRICATOR® (IdeS)
Manufacturer	Genovis AB, Lund, Sweden
Appearance	Colourless, clear solution
Buffer	10 mM sodium phosphate, 140mM NaCl, pH 7.4 with no preservatives added.
Storage Temperature	-20°C. Avoid freeze/thaw cycles
Shipping Conditions	Shipped cold with gel ice refrigerants
Quality Control	Detect 10 ng FabRICATOR at dilution 1:500 using Western blot chromogenic detection