

Product Name	FabRICATOR® protease (IdeS), lyophilized
Product Number	A0-FR1-008
Product composition	800 (8×100) Units* of FabRICATOR® proteolytic enzyme lyophilized from phosphate buffer saline (PBS).
Physical purity	Purified to ≥ 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.
Manufacturer	Genovis AB, Sweden
Storage temp	-20°C
Appearance	White to light yellow powder
Source	Cloned from <i>Streptococcus pyogenes</i> and overexpressed in <i>E. coli</i> .
pH in PBS	6.5-7.5
Sterility	Sterile according to tests**
Stability	2 years in -20°C

\*Unit Definition      One unit is defined as the amount of enzyme required to fragment ≥95% of 1 µg of human IgG in 10 mM sodium phosphate, 137 mM NaCl, 2.7mM KCl pH 7.4 at 37°C for 30 min as monitored by SDS-PAGE.

\*\*Sterile tests      Tested on Sabouraud Dextrose Agar, Fluid Thioglycollate Medium and Blood Agar.

---

Telephone:	+46 (0)46 10 12 30	Visiting address:	Postal address:	VAT number:	SE556574534501
Telefax:	+46 (0)46 12 80 20	Scheelevägen 8	Box 790		
E-mail:	info@genovis.com	Medicon Village	SE-223 07 Lund		
Internet:	www.genovis.com	SE-223 63 Lund	Sweden		
		Sweden	Sweden		