

Anti-FabRICATOR[®] antibody

Product no
A3-AF1-010 100 µl


Properties

Product Name	Anti-FabRICATOR [®] antibody
Description	Goat polyclonal antibody to FabRICATOR [®]
Purity	Protein G purified
Clonality	Polyclonal
Isotype	IgG
Concentration	4 mg/ml
Immunogen	FabRICATOR [®]
Appearance	Form: Liquid Clarity: Clear Color: Colorless Odor: None
Buffer	Anti-FabRICATOR [®] antibody is supplied in solution of 10 mM sodium phosphate, 150 mM NaCl, pH 7.4 with no preservative added.
Storage	Anti-FabRICATOR [®] antibody is shipped on 4°C / 39°F. Upon delivery aliquot and store at -20°C / -4°F. Avoid freeze/thaw cycles.

Applications

Western blot

The antibody is tested for use in western blot with chromogenic detection. Use a dilution of 1:1000 to 1:500.



Anti-FabRICATOR[®] 1:500 dilution

FabRICATOR[®] 10 ng (0.5 units)

Predicted band size: 37 kDa

Anti-FabRICATOR[®] antibody is for R&D use only.