

Product Name	GingiskHAN™ (KGP), lyophilized
Product Number	B0-GKH-020
Product composition	2,000 Units* of GingiskHAN™ proteolytic enzyme lyophilized from Tris buffer.
Physical purity	Purified to $\geq 80\%$ homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.
Manufacturer	Genovis AB, Sweden
Storage temp	-20°C
Appearance	White to light yellow powder
Source	Purified from <i>Porphyromonas gingivalis</i> .
pH in Tris	8.0
Sterility	Sterile according to tests**
Stability	1 year in -20°C

*Unit Definition One unit cleaves $\geq 95\%$ of 1 μ g human IgG1 when incubated in 0.1 M Tris, pH 8.0 at 37°C for 1 hour.

**Sterile tests Tested on Sabouraud Dextrose Agar, Fluid Thioglycollate Medium and Blood Agar.

Telephone:	+46 (0)46 10 12 30	Visiting address:	Postal address:	VAT number:	SE556574534501
Telefax:	+46 (0)46 12 80 20	Scheelevägen 8	Box 790		
E-mail:	info@genovis.com	Medicon Village	SE-223 07 Lund		
Internet:	www.genovis.com	SE-223 63 Lund	Sweden		
		Sweden	Sweden		