

## Product Specification

Version 01  
Issued June 2026

<b>Product Name</b>	AlligBAITOR™ Lyophilized
<b>Product Number</b>	A4-AB1-020
<b>Product Composition</b>	The AlligBAITOR is lyophilized in TBS, pH 7.6, with a total of 2000 units per vial.
<b>Manufacturer</b>	Genovis AB, Sweden
<b>Appearance</b>	White to light yellow powder
<b>Binding capacity</b>	≥90% complex formation between AlligBAITOR and IgG (1U/ug) when incubated in PBS, pH 7.4, at room temperature for 15 minutes with no detectable free AlligBAITOR, based on SEC analysis.
<b>Physical Purity</b>	≥90% homogeneity based on SDS-PAGE analysis using Coomassie Blue detection.
<b>Source</b>	AlligBAITOR is recombinantly expressed in <i>E. coli</i> .
<b>Bioburden*</b>	Free from microbial contamination
<b>Quality Standard</b>	ISO 9001:2015
<b>Shipping Conditions</b>	Shipped cold with gel ice refrigerants
<b>Storage Temperature</b>	-20°C

\* Tested on Dextrose Agar, Fluid Thioglycollate Medium and Blood Agar.