

Product Specification

Version 02
Issued April 2026

Product Name	FabDELLO® Lyophilized
Product Number	B1-BD1-008, B1-BD1-096
Product Composition	<ul style="list-style-type: none">• The FabDELLO enzyme (B0-BD1-008/B0-BD1-096) is lyophilized in TBS with a total of 100 units per plate vial.• 200 mM CaCl₂ solution in water (CACL2-102)
Manufacturer	Genovis AB, Sweden
Enzyme Appearance	White to light yellow powder
Enzyme Activity	One unit FabDELLO digests ≥90% of 1 µg human IgG1 when incubated in TBS with 10 mM CaCl ₂ , pH 7.6, at 37°C, for 2 hours based on SDS-PAGE analysis.
Purity	Absence of contaminating activity as determined by SDS-PAGE analysis using Coomassie Blue detection.
Source	The FabDELLO enzyme is recombinantly expressed in <i>E. coli</i> .
Quality Standard	ISO 9001:2015
Shipping Conditions	Shipped cold with gel ice refrigerants
Storage Temperature	The FabDELLO enzyme is stored at -20°C. The CaCl ₂ solution is stored at +4-8 °C.