



Last revised Jan 2023

<b>Product Name</b>	FabRICATOR® HPLC
<b>Product Number</b>	A0-FRC-050
<b>Product Composition</b>	FabRICATOR enzyme covalently immobilized on HPLC compatible resin packed in peek column 2.1 x 50mm in PBS, pH 7.4 with 0.02% sodium azide.
<b>Manufacturer</b>	Genovis AB, Sweden
<b>Source</b>	FabRICATOR (IdeS) enzyme is cloned from <i>Streptococcus pyogenes</i> and expressed in <i>E. coli</i> .
<b>Column performance</b>	≥ 95% digestion of 10 µg hIgG1 in less than 5 minutes residence time at 37°C as monitored by RP-LC.
<b>Shipping Conditions</b>	Product shipped cold
<b>Storage Temperature</b>	+4-8°C