



Issued March 2021

<b>Product Name</b>	FucosEXO™
<b>Product Number</b>	G1-FM1-020
<b>Product Composition</b>	2000 Units of FucosEXO, exoglycosidase lyophilized in Tris buffer saline (TBS).
<b>Manufacturer</b>	Genovis AB, Sweden.
<b>Appearance</b>	White to light yellow powder.
<b>Source</b>	FucosEXO enzymes are recombinantly expressed in <i>E.coli</i> .
<b>Bioburden</b>	Absence of microbial contamination. Tested on Sabouraud Dextrose Agar, Fluid Thioglycollate Medium and Blood Agar.
<b>Enzyme Activity</b>	One unit of FucosEXO hydrolyzes $\alpha$ 1-2 and $\alpha$ 1-3-linked fucose from $\geq 90\%$ of 11 nmol 2'-fucosyllactose and 4.5 nmol 3-fucosyllactose when incubated in 20 mM Tris pH 6.8 at 37°C for 30 minutes.
<b>Physical Purity</b>	Purified to $\geq 95\%$ homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.
<b>Shipping Conditions</b>	Shipped with gel ice refrigerants.
<b>Storage Temperature</b>	-20°C