

Product Specification

Version 02
Issued December 2025

Product Name	GalNAcEXO™ Lyophilized
Product Number	G1-NA1-020
Product Composition	The GalNAcEXO enzyme is lyophilized in TBS, pH 7.6, with a total of 2000 units per vial.
Manufacturer	Genovis AB, Sweden
Appearance	White to light yellow powder
Enzyme Activity	One unit of GalNAcEXO hydrolyzes $\geq 95\%$ α -linked GalNAc residues of 2 nmol 4-Nitrophenyl 2-acetamido-2-deoxy- α -D-galactopyranoside, when incubated in 20 mM Tris, 1% ethanol, pH 6.8, at 37°C, for 30 minutes based on RP-LC analysis.
Physical Purity	Purified to $\geq 95\%$ homogeneity based on SDS-PAGE analysis with Coomassie Blue detection.
Source	The GalNAcEXO enzymes are recombinantly expressed in <i>E. coli</i> .
Bioburden*	Free from microbial contamination
Quality Standard	ISO 9001:2015
Shipping Conditions	Shipped cold with gel ice refrigerants
Storage Temperature	-20°C

* Tested on Dextrose Agar, Fluid Thioglycollate Medium and Blood Agar.