

Product Specification

Version 01
Issued March 2026

Product Name	SEQguard™ Dino Preserve
Product Number	PS-SDP-024
Product Composition	SEQguard Dino Preserve is supplied as an aqueous solution free from any toxic substances. 150 µl solution per vial sufficient to stabilize 24 purified RNA samples when used according to Plasmidsaurus' usage guidelines.
Supplier	Genovis AB, Sweden
Appearance	Clear colorless solution
Usage	Transport medium for ambient temperature shipment of purified RNA samples.
Quality Standard	ISO 9001:2015
Shipping Conditions	Shipped at ambient temperature
Storage Temperature	15-30°C. It also tolerates freeze-thaw cycles.